

SOX10^(M)


The Cleanest SOX




The Cleanest SOX

SOX10 (M) [BC34] from Biocare Medical is an exclusive IVD mouse monoclonal antibody that provides a highly specific and sensitive antibody for immunohistochemistry and is virtually background free.

SOX10 has been shown to stain the vast majority of cutaneous and metastatic melanomas, including 97-100% of spindle cell and desmoplastic melanomas, 100% of nevi, the majority of oligodendrogliomas and a large percentage of astrocytomas and poorly differentiated glioblastomas. PATENT PENDING.



Melanoma stained with SOX10 (M)



Spindle Cell Melanoma stained with SOX10 (M)

Ordering Information	Concentrate	Predilute	VP
SOX10 (M)	ACI 3099 A, C	API 3099 AA	AVI 3099 G

www.biocare.net

1. Mohamed A, et al. Appl Immunohistochem Mol Morphol. 2013 Dec; 21(6):506-10. 2. Pusch C, et al. Hum Genet. 1998 Aug; 103(2):115-23. 3. Mollaaghababa R, Pavan WJ. Oncogene. 2003 May; 22(20):3024-34. 4. Bondurand N, et al. FEBS Lett. 1998 Aug; 432(3):168-72. 5. Bannaykh SI, et al. J Neurooncol. 2006 Jan; 76(2):115-27. 6. Britsch S, et al. Genes Dev. 2001 Jan; 15(1):66-78. 7. Feng Z, et al. J Cutan Pathol. 2011 Aug; 38(8):616-24.